

Certificate of Analysis

Product name: Anti-T (Arachis)
Blood Grouping Reagent

Product code(s): 03-009 – 03-014, 03-029 – 03-037

Isotype: Lectin Extract

Batch Numbers: 14150402

Date Manufactured: 30/04/2015

Expiry Date: 29/04/2017

Appearance: Clear – Yellow Coloured liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	6.4 – 7.2	7.03
Conductivity (µs/cm)	N/A	N/A
Osmolarity (mOsm/kg)	290-1000	484

Titre: This is a Slide or tube reacting reagent.

Avidity: Matches at least the minimum requirement of 1:4
(Spin/Tube at 23°C).

Preservatives used: Sodium Azide (0.09% [w/v])

Sterility: NA

Storage: +2°C to +8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Sodium azide is poisonous and is an irritant. It can be fatal if swallowed and may cause cardiac disturbances. Although present in these reagents in small quantities, appropriate health and safety interventions are recommended. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Quality Assurance:

Date: 07 May 2015

