

Certificate of Analysis

Product name: Filtered Reagent Grade Water

Product code(s): 04-036

Batch Numbers: 18160804

Date Manufactured: 23 Aug 2016

Expiry Date: 22 Aug 2017

Appearance: Clear liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	<10	6.50
Osmolarity(mOsm/kg)	N/A	N/A

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Store at room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no health and safety issues associated with this product.

Signed for Quality Assurance:

Date: 08 Sep 2016

