

## Certificate of Analysis

**Product name: Fixol**  
Fixative reagent for the FCD kit

**Product code(s): 05-051**

**Batch Numbers: 51160505**

**Date Manufactured: 06 May 2016**

**Expiry Date: 05 May 2018**

**Appearance: Clear Fluid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	8.50-9.50	8.72
<b>Conductivity(<math>\mu</math>s/cm)</b>	1750-2150	1846
<b>Osmolarity(mOsm/kg)</b>	N/A	N/A

**Preservatives used:** Alcohol Based

**Sterility:** Not required

**Storage:** Room temperature

**Country of Origin:** United Kingdom

### **Health and Safety:**

#### *Precautionary Notes*

Methanol is a very toxic substance if swallowed. It can be fatal if swallowed. It is also toxic by inhalation. It is an irritant to the eyes.

Methanol is highly flammable.

It is strongly recommended that the user holds the appropriate COSHH or MSDS information on this product.

**Signed for Quality Assurance:**

**Date: 10 May 2016**

