

Certificate of Analysis

Product name: Antibiotic Concentrate

Product code(s): 07-011 and 07-012

Batch Numbers: 82170103

Date Manufactured: 31 Jan 17

Expiry Date: 30 Jan 18

Appearance: Yellow coloured liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	3.5-5.5	4.39
Conductivity(μ s/cm)	N/A	N/A
Osmolarity(mOsm/kg)	N/A	N/A

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade B clean room and tested for sterility.

Storage: Store Frozen

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Gentamicin sulphate is classified as 'Harmful'. While this antibiotic is dissolved in the reagent, appropriate health and safety interventions are recommended. Exposure of preserved reagents directly to the skin and eyes should be avoided. Appropriate Personal Protective Equipment (PPE) should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Release:

Date: 10 Feb 2017

