

## Certificate of Analysis

**Product name: Cryoprotectant Reagent (CPR)**

**Product code(s): 04-211 – 04-214, 04-165, and 04-210PL**

**Batch Numbers: 61160705**

**Date Manufactured: 15 Jul 16**

**Expiry Date: 14 Jul 18**

**Appearance: Clear Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	3.50-5.50	4.04
<b>Conductivity(us/cm)</b>	4200-4850	4350
<b>Osmolarity(mOsm/kg)</b>	880-970	911

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility.

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

There are no specific Health and Safety problems associated with the use of CPR.

**Signed for Quality Assurance:**

**Date: 20 Jul 2016**

