

Certificate of Analysis

Product name: Cryoprotectant Reagent (CPR)

Product code(s): 04-211 – 04-214, 04-165, and 04-210PL

Batch Numbers: 61161101

Date Manufactured: 09 Nov 16

Expiry Date: 08 Nov 18

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	3.50-5.50	3.51
Conductivity(us/cm)	4200-4850	4553
Osmolarity(mOsm/kg)	880-970	881

Preservatives used: None

Sterility: Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of CPR.

Signed for Quality Assurance:

Date: 21 Nov 2016

