

Certificate of Analysis

Product name: DCD Cryo Medium
(Decellular Dermis Cryo Medium)

Product code(s): 04-351

Batch Numbers: 83150403

Date Manufactured: 14/04/2015

Expiry Date: 13/04/2016

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	7.10-7.50	7.25
Conductivity(μs/cm)	790-850	822
Osmolarity(mOsm/kg)	2400-2520	2520

Meets the requirements set out by Tissue Services, Liverpool.

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade C clean room and tested for sterility.

Storage: Can be stored at room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of DCD Cryo Medium.

Signed for Release:

Date: 22 Apr 2015

