

## Certificate of Analysis

**Product name: Frozen Cell Recovery Solution (FCRS)**

**Product code(s): 04-161 – 04-162**

**Batch Numbers: 41160903**

**Date Manufactured: 27 Sep 2016**

**Expiry Date: 26 Sep 2018**

**Appearance: Clear Liquid**

**Quality Control:**

	Acceptable Ranges	Batch result
<b>pH</b>	6.40-6.80	6.40
<b>Conductivity(<math>\mu</math>s/cm)</b>	44000-52000	46000
<b>Osmolarity(mOsm/kg)</b>	1080-1120	1120

**Preservatives used:** EDTA (anti-fungal)

**Sterility:** Terminally filtered with 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

This reagent contains a high concentration of sodium chloride but presents only a very minor risk to staff handling the solution.

**Signed for Quality Assurance:**

**Date: 10 Oct 2016**

