

## Certificate of Analysis

### Product name: Frozen Cell Recovery Solution (FCRS)

Product code(s): 04-161 – 04-162

Batch Numbers: 41161104

Date Manufactured: 15 Nov 2016

Expiry Date: 14 Nov 2018

Appearance: Clear Liquid

#### Quality Control:

	Acceptable Ranges	Batch result
pH	6.40-6.80	6.40
Conductivity( $\mu$ s/cm)	44000-52000	51033
Osmolarity(mOsm/kg)	1080-1120	1090

Preservatives used: EDTA (anti-fungal)

Sterility: Terminally filtered with 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

Storage: 2-8°C

Country of Origin: United Kingdom

#### Health and Safety:

##### *Precautionary Notes*

This reagent contains a high concentration of sodium chloride but presents only a very minor risk to staff handling the solution.

Signed for Quality Assurance:

Date: 29 Nov 2016

