

Certificate of Analysis

Product name: Filtered Reagent Grade Water

Product code(s): 04-036

Batch Numbers: 18161203

Date Manufactured: 08 Dec 2016

Expiry Date: 07 Dec 2017

Appearance: Clear liquid

Quality Control:

| | Acceptable Ranges | Batch result |
|---------------------------|-------------------|--------------|
| pH | N/A | N/A |
| Conductivity(μ s/cm) | <10 | 7.14 |
| Osmolarity(mOsm/kg) | N/A | N/A |

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade B clean room and tested for sterility.

Storage: Store at room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no health and safety issues associated with this product.

Signed for Quality Assurance:

Date: 23 Dec 2016

