

## Certificate of Analysis

**Product name: Fixol**  
Fixative reagent for the FCD kit

**Product code(s): 05-051**

**Batch Numbers: 51160916**

**Date Manufactured: 07 Sep 2016**

**Expiry Date: 30 Aug 2018**

**Appearance: Clear Fluid**

**Quality Control:**

|   | <b>Acceptable Ranges</b> | <b>Batch result</b> |
|---|--------------------------|---------------------|
| <b>pH</b>                                 | 8.50-9.50                | 8.55                |
| <b>Conductivity(<math>\mu</math>s/cm)</b> | 1750-2150                | 2050                |
| <b>Osmolarity(mOsm/kg)</b>                | N/A                      | N/A                 |

**Preservatives used:** Alcohol Based

**Sterility:** Not required

**Storage:** Room temperature

**Country of Origin:** United Kingdom

### **Health and Safety:**

#### *Precautionary Notes*

Methanol is a very toxic substance if swallowed. It can be fatal if swallowed. It is also toxic by inhalation. It is an irritant to the eyes.

Methanol is highly flammable.

It is strongly recommended that the user holds the appropriate COSHH or MSDS information on this product.

**Signed for Quality Assurance:**

**Date: 26 Sep 2016**

