

Certificate of Analysis

Product name: Heart Valve Cryomedia Solution (HVCS) SNBTS

Product code(s): 04-321

Batch Numbers: 19160602

Date Manufactured: 28 Jun 2016

Expiry Date: 27 Jun 2017

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	This batch result
pH	7.20-7.40	7.38
Conductivity(μs/cm)	9000-11000	9205
Osmolarity(mOsm/L)	1950-2150	1984

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in grade B clean room and product tested for sterility using various microbiological methods before clearance.

Storage: Frozen

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Dimethyl Sulfoxide (DMSO) is classed as an 'Irritant' and care should be taken with this reagent as it contains 15% DMSO. It is readily absorbed through the skin and skin contamination should be washed off immediately with copious amounts of water.

Signed for Quality Assurance:

Date: 03 Jun 2016

