

Certificate of Analysis

Product name: Hypotonic Buffer Base

Product code(s): 04-333

Batch Numbers: 76160704

Date Manufactured: 19 Jul 2016

Expiry Date: 18 Jul 2017

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	7.80-8.20	7.89
Conductivity(μs/cm)	700-1100	1083
Osmolarity(mOsm/kg)	20-30	25

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade B clean room and tested for sterility.

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Both EDTA and TRIS are 'Irritants' and every care should be taken when handling the reagent. Although present in very small quantities, appropriate health and safety interventions are recommended. Exposure of the reagent directly to the skin and eyes should be avoided. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Quality Assurance:

Date: 25 Jul 2016

