

Certificate of Analysis

Product name: Isotonic Saline

Product code(s): 04-025 – 04-027

Batch Numbers: 66160902

Date Manufactured: 07 Sep 2016

Expiry Date: 06 Sep 2017

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	N/A	N/A
Conductivity(μs/cm)	N/A	N/A
Osmolarity(mOsm/kg)	285-295	289

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade B clean room and tested for sterility.

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of Isotonic Saline.

Signed for Quality Assurance:

Date: 26 Sep 2016

