

Certificate of Analysis

Product name: LISSC (Low Ionic Strength Saline Concentrate)

Product code(s): 04-019 – 04-021, 04-122, and 04-122

Batch Numbers: 39160802

Date Manufactured: 31 Aug 2016

Expiry Date: 30 Aug 2018

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	6.50-6.90	6.80
Conductivity(μ s/cm)	3400-4000	3820
Osmolarity(mOsm/kg)	285-295	293

Preservatives used: Sodium Azide (0.09% [w/v])

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Sodium azide is poisonous and is an irritant. It is fatal if swallowed and may cause cardiac disturbances. Although present in these reagents in very small quantities, appropriate health and safety interventions are recommended. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Quality Assurance:

Date: 12 Sep 2016

