

Certificate of Analysis

Product name: PBSS pH7.00
(Phosphate Buffered Serology Saline Azide Free)

Product code(s): 04-101 – 04-104

Batch Numbers: 06161205

Date Manufactured: 14 Dec 2016

Expiry Date: 13 Dec 2018

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	This batch result
pH	6.95-7.05	6.96
Conductivity(μs/cm)	N/A	N/A
Osmolarity(Osm/kg)	285-295	290

Preservatives used: None

Sterility: Terminally filtered with 0.2 μ m filter in a grade B clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of PBSS pH7.00.

Signed for Quality Assurance:

Date: 23 Dec 2015

