

Certificate of Analysis

Product name: SkinCryo Medium

Product code(s): 04-336

Batch Number: 73160602

Date Manufactured: 30 Jun 16

Expiry Date: 29 Jun 17

Appearance: Bright Red Liquid

Quality Control:

| | Acceptable Ranges | Batch result |
|----------------------------|--------------------------|---------------------|
| pH | 7.10-7.50 | 7.45 |
| Conductivity(ms/cm) | 7.90-8.50 | 7.95 |
| Osmolarity(mOsm/kg) | 2200-2550 | 2475 |

Preservatives used: None

Sterility: Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of SkinCryo Medium

Signed for Quality Assurance:

Date: 11 Jul 2016

