

## Certificate of Analysis

**Product name: Red Cell Support Solution (RCSS)**

**Product code(s): 04-131-04-132 and 04-135-04-138**

**Batch Numbers: 35160806**

**Date Manufactured: 23 Aug 2016**

**Expiry Date: 22 Aug 2018**

**Appearance: Clear Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	6.80-7.20	6.88
<b>Conductivity(<math>\mu</math>s/cm)</b>	3100-3300	3255
<b>Osmolarity(mOsm/kg)</b>	300-315	311

**Preservatives used:** Antibiotics: Sulfamethoxazole and Trimethoprim

**Sterility:** Terminally filtered with 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

There are no specific Health and Safety problems associated with the use of RCSS. However, users should apply GLP when using any reagent.

**Signed for Quality Assurance:**

**Date: 07 Sep 2016**

