

Certificate of Analysis

Product name: Red Cell Support Solution (RCSS)

Product code(s): 04-131-04-132 and 04-135-04-138

Batch Numbers: 35170107

Date Manufactured: 01 Feb 2017

Expiry Date: 31 Jan 2019

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	6.80-7.20	6.87
Conductivity($\mu\text{s}/\text{cm}$)	3100-3300	3300
Osmolarity(mOsm/kg)	300-315	309

Preservatives used: Antibiotics: Sulfamethoxazole and Trimethoprim

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8 $^{\circ}\text{C}$

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of RCSS. However, users should apply GLP when using any reagent.

Signed for Quality Assurance:

Date: 15 Feb 2017

