

## Certificate of Analysis

**Product name:** Reagent Grade Water

**Product code(s):** 04-035

**Batch Numbers:** 17170362

**Date Manufactured:** 10 Mar 17

**Expiry Date:** 09 Mar 18

**Appearance:** Clear

**Quality Control:**

	Acceptable Ranges	Batch result
pH	N/A	N/A
Conductivity( $\mu\text{s}/\text{cm}$ )	<10	9.63
Osmolarity(mOsm/kg)	N/A	N/A

**Preservatives used:** None

**Sterility:** No chemicals used in this preparation. Water is run through the reverse osmosis unit and collected without it being stored. It is collected straight into the appropriate bags or bottles with as little exposure to the atmosphere as possible.

**Storage:** Store at room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**  
*Precautionary Notes*

**Signed for Release:**

**Date:** 10 Mar 2017

