

## Certificate of Analysis

**Product name: SkinCryo Medium**

**Product code(s): 04-336**

**Batch Number: 73150705**

**Date Manufactured: 28/07/15**

**Expiry Date: 27/07/16**

**Appearance: Bright Red Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	7.10-7.50	7.30
<b>Conductivity(ms/cm)</b>	7.90-8.50	7.95
<b>Osmolarity(mOsm/kg)</b>	2200-2550	2301

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2µm filter in a grade C clean room and tested for sterility.

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

There are no specific Health and Safety problems associated with the use of SkinCryo Medium

**Signed for Quality Assurance:**

**Date: 10 Aug 2015**

