

Certificate of Analysis

Product name: TRIS Buffered Isotonic Saline(TBIS)

Product code(s): 04-141 – 04-144

Batch Numbers: 37170303

Date Manufactured: 08 Mar 2017

Expiry Date: 07 Mar 2018

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	7.30-7.60	7.56
Conductivity(ms/cm)	14-16.5	16.01
Osmolarity(mOsm/kg)	285-300	294

Preservatives used: None

Sterility: Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility.

Storage: Can be stored at room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific health and safety problems associated with the use of TBIS.

Signed for Quality Assurance:

Date: 15 Mar 2017

