

## Certificate of Analysis

**Product name: Wash Buffer Base**

**Product code(s): 04-331**

**Batch Numbers: 76160802**

**Date Manufactured: 09 Aug 2016**

**Expiry Date: 08 Aug 2017**

**Appearance: Clear Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	7.80-8.20	7.96
<b>Conductivity(<math>\mu</math>s/cm)</b>	1510-1600	1569
<b>Osmolarity(mOsm/kg)</b>	290-310	297

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2 $\mu$ m filter in a grade C clean room and tested for sterility.

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

### **Health and Safety:**

#### *Precautionary Notes*

EDTA is considered an 'Irritant' and every care should be taken when handling the reagent. Although present in very small quantities, appropriate health and safety interventions are recommended. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

**Signed for Quality Assurance:**

**Date: 16 Aug 2016**

