

Certificate of Analysis

Product name: Cryoprotectant Reagent (CPR)

Product code(s): 04-211 – 04-214, 04-165, and 04-210PL

Batch Numbers: 61170406

Date Manufactured: 12 Apr 17

Expiry Date: 11 Apr 19

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	5.80-6.70	6.13
Conductivity(us/cm)	4200-4850	4550
Osmolarity(mOsm/kg)	880-970 (1:5*)	927

*Dilution made 1 part water, 4 parts Cryoprotectant Reagent.

Preservatives used: None

Sterility: Terminally filtered with 0.2µm filter in a grade B clean room and tested for sterility.

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

There are no specific Health and Safety problems associated with the use of CPR.

Signed for Quality Assurance:

Date:

26 Apr 2017

