

Certificate of Analysis

Product name: Checkzyme (*Glycine soja* extract)

Product code(s): 04-037, 04-037/100, 04-038 and 04-151

Batch Numbers: 15170405

Date Manufactured: 13 Apr 2017

Expiry Date: 12 Apr 2019

Appearance: Clear- Yellow coloured liquid

Quality Control:

	Acceptable Ranges	This batch result
pH	6.9-7.1	7.10
Conductivity(μ s/cm)	N/A	N/A
Osmolarity(mOsm/kg)	285-295	287

Preservatives used: Sodium azide (0.09% [w/v])

Storage: 2-8°C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Care should be taken if glass vials are present in this kit. In the case of breakage resulting in a sharps injury, the appropriate first aid measures should be taken.

Sodium azide is poisonous and is an irritant. It can be fatal if swallowed and may cause cardiac disturbances. Although present in these reagents in small quantities, appropriate health and safety interventions are recommended. Exposure of preserved reagents directly to the skin and eyes should be avoided. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same time.

Signed for Release:

Date: 26 Apr 2017

