

## Certificate of Analysis

**Product name: SkinCryo Medium**

**Product code(s): 04-336**

**Batch Number: 73170402**

**Date Manufactured: 27 Apr 17**

**Expiry Date: 26 Apr 18**

**Appearance: Red Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	7.10-7.50	7.22
<b>Conductivity(<math>\mu</math>s/cm)</b>	7900-8500	8443
<b>Osmolarity(mOsm/kg)</b>	2200-2550	2390

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

There are no specific Health and Safety problems associated with the use of SkinCryo Medium

**Signed for Quality Assurance:**

**Date: 04 May 2017**

