

## Certificate of Analysis

**Product name: Frozen Cell Recovery Solution (FCRS)**

**Product code(s): 04-161 BAG & 04-162 BAG**

**Batch Numbers: 41170508**

**Date Manufactured: 23 May 2017**

**Expiry Date: 22 May 2019**

**Appearance: Clear Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	6.40-6.80	6.52
<b>Conductivity(<math>\mu</math>s/cm)</b>	44000-52000	51567
<b>Osmolarity(mOsm/kg)</b>	1080-1120	1095

**Preservatives used:** EDTA (anti-fungal)

**Sterility:** Terminally filtered with a 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

### **Health and Safety:**

#### *Precautionary Notes*

This reagent contains a high concentration of sodium chloride but presents only a very minor risk to staff handling the solution.

Find the Safety Data sheet on our website at [www.sourcebioscience.com](http://www.sourcebioscience.com)

**Signed for Quality Assurance:**

**Date: 05 Jun 2017**

