

Certificate of Analysis

Product name: Fetal Cell Detection Kit (FCD)

Product code(s): 05-050 – 05-052

Batch Numbers: 05170522 A/B

Date Manufactured: Eluting: 17 May 2017 / Counterstain: 22 May 2017

Kit Expiry Date: 16 May 2019

Appearance: Eluting Reagent Black/Deep Purple Fluid
Counterstain Reagent Bright Red Fluid

Quality Control:	pH	Range	This Batch
	Eluting Reagent:		1.80-2.10
Counterstain Reagent:		N/A	N/A

The performance of the above components has been assessed by Kleihauer testing of negative and positive control samples. All tests met the performance specifications for this product and passed the assessment.

[Demonstration of fetal hemoglobin in erythrocytes of a blood smear].

KLEIHAUER E, BRAUN H, BETKE K
Klin Wochenschr 1957 (35):637-638

Preservatives used: Eluting Reagent Alcohol Based
Counterstain Reagent 0.10% w/v Sodium azide

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety: *Precautionary Notes*

Find the Safety Data sheet on our website at www.sourcebioscience.com

Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests.
Seek medical advice if any of these reagents are swallowed.

Signed for Quality Assurance:

Date: 09 Jun 2017

