

## Certificate of Analysis

**Product name: Red Cell Support Solution (RCSS)**

**Product code(s): 04-131 BAG, 04-132 BAG**

**Batch Numbers: 35170609**

**Date Manufactured: 29 Jun 2017**

**Expiry Date: 28 Jun 2019**

**Appearance: Clear Liquid**

**Quality Control:**

	<b>Acceptable Ranges</b>	<b>Batch result</b>
<b>pH</b>	6.80-7.20	6.94
<b>Conductivity(<math>\mu</math>s/cm)</b>	3100-3300	3260
<b>Osmolarity(mOsm/kg)</b>	300-315	304

**Preservatives used:** Antibiotics: Sulfamethoxazole and Trimethoprim

**Sterility:** Terminally filtered with 0.2 $\mu$ m filter in a grade B clean room and tested for sterility.

**Storage:** 2-8°C

**Country of Origin:** United Kingdom

**Health and Safety:**

*Precautionary Notes*

Find the Safety Data sheet on our website at [www.sourcebioscience.com](http://www.sourcebioscience.com)

Use appropriate Personal Protective Equipment when using this reagent.

**Signed for Quality Control & Batch Release:**

**Date: 10 July 2017**

