

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Isotonic Saline

Product code(s): 04-025

Batch Number(s): 66170803

Date Manufactured: 14 Aug 17

Expiry Date: 13 Aug 18

Appearance: Clear liquid

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	289

Preservatives used: None

Sterility: Terminally filtered through a 0.2 μm filter in a grade B clean room and sterility tested.

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Precautionary Notes

There are no specific Health and Safety problems associated with the use of Isotonic Saline.

Signed for Quality Control: 

Date: 23 Aug 17

Signed for Certificate Check: 

Date: 23 Aug 17