

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** PBSS AF ( Phosphate Buffered Serology Saline, Azide Free) pH 7.30

**Product code(s):** 04-046 and 04-047

**Batch Number(s):** 25170705

**Date Manufactured:** 01 Aug 2017

**Expiry Date:** 31 Jul 2019

**Appearance:** Clear Liquid

**Quality Control:**

	Specification	This batch result
pH	7.25 - 7.35	7.34
Conductivity( $\mu\text{s}/\text{cm}$ )	N/A	N/A
Osmolality(mOsm/kg)	285 - 295	295

**Preservatives used:** None

**Sterility:** Terminally filtered with a 0.2 $\mu\text{m}$  sterile filter in a grade B clean room and tested for sterility.

**Storage:** Room Temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

09 Aug 17

**Signed for Certificate Check:**



**Date:**

16 Aug 17