

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Heart Valve Cryomedia Solution (HVCS)

**Product code(s):** 04-311

**Batch Number(s):** 71170705

**Date Manufactured:** 25 Jul 2017

**Expiry Date:** 24 Jul 2018

**Appearance:** Red Translucent liquid

### Quality Control:

	Specification	This batch result
pH	7.00 - 7.50	7.33
Conductivity( $\mu\text{s}/\text{cm}$ )	8500 - 12000	9233
Osmolality(mOsm/kg)	2600 - 3000	2799

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in grade B clean room and tested for sterility

**Storage:** Room Temperature

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

It is advised that Personal Protective Equipment (PPE) be worn when handling this reagent.  
Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Signed for Certificate Check:**



**Date:**

09 Aug 17

**Date:**

09 Aug 17