

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsever's Solution  
**Product code(s):** 04-182  
**Batch Number(s):** 32170712  
**Date Manufactured:** 03 Jul 2017  
**Expiry Date:** 02 Jul 2018  
**Appearance:** Clear liquid

### Quality Control:

	Specification	This batch result
pH	6.80 - 7.30	7.06
Conductivity( $\mu\text{s}/\text{cm}$ )	11500 - 13000	12433
Osmolality(mOsm/kg)	335 - 360	343

**Preservatives used:** Chloramphenicol and Neomycin sulphate

**Sterility:** Terminally filtered through a 0.2 $\mu\text{m}$  filter in a grade B clean room and sterility tested.

**Storage:** 2 - 8 °C

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Appropriate Personal Protective Equipment (PPE) should be worn at all times whilst carrying out tests. Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

12 Jul 2017

**Signed for Certificate Check:**



**Date:**

13 Jul 17