

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** LISS (Low Ionic Strength Saline)  
**Product code(s):** 04-016 TRIAL  
**Batch Number(s):** 03170816  
**Date Manufactured:** 17 Aug 2017  
**Expiry Date:** 16 Aug 2018  
**Appearance:** Clear liquid

**Quality Control:**

|   | Specification | This batch result |
|---|---------------|-------------------|
| pH                                      | 6.50 - 6.90   | 6.74              |
| Conductivity( $\mu\text{s}/\text{cm}$ ) | 3500 - 3800   | 3610              |
| Osmolality(mOsm/kg)                     | 285 - 295     | 294               |

**Preservatives used:** Sodium azide (0.09% [w/v])

**Sterility:** N/A

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Wear appropriate Personal Protective Equipment (PPE) when handling this reagent.  
Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

17 Aug 17

**Signed for Certificate Check:**



**Date:**

21 Aug 17