

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Reagent Grade Water

**Product code(s):** 04-035

**Batch Number(s):** 17170871

**Date Manufactured:** 10 Aug 2017

**Expiry Date:** 09 Aug 2018

**Appearance:** Clear liquid

**Quality Control:**

|   | Specification | This batch result |
|---|---------------|-------------------|
| pH                                      | N/A           | N/A               |
| Conductivity( $\mu\text{s}/\text{cm}$ ) | 0 - 10        | 3.24              |
| Osmolality(mOsm/kg)                     | N/A           | N/A               |

**Preservatives used:** No preservatives used

**Sterility:** Water is UV treated and processed into gamma irradiated bags. Although exposure to the environment is limited, we do not guarantee the sterility of the water.

**Storage:** Store at Room Temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

There are no H&S implications from this product;

Please download the safety data sheet from [www.sourcebioscience.com](http://www.sourcebioscience.com)

**Signed for Quality Control:**



**Date:**

10 AUG 17

**Signed for Certificate Check:**



**Date:**

11 Aug 17