

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: POLYAGGLUTINATION KIT
Product code(s): 03-061
Batch Number(s): 79161006
Date Manufactured: 19 Oct 16
Expiry Date: 18 Aug 18
Appearance: Clear – pale brown liquid
Kit made up of four lectin extracts: 1. Arachis hypogea - 14160903, 2. Glycine soja - 15160803
3. Salvia Sclarea - 49160905, 4. Dolichos biflorus - 29161010

Quality Control:

	Specification	This batch result
pH	N/A	N/A
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	N/A	N/A

No specific guidelines are laid down for a polyagglutination kit but the kit is expected to show a pattern of agglutination with affected human red cells which will classify them as being T, Tn, Tk or Cad. It may be possible to further determine Th and Tx affected red cells

Preservatives used: Sodium azide (0.1% [w/v])

Sterility: N/A

Storage: 2 - 8 °C

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Sodium azide is poisonous and is an irritant. It is fatal if swallowed and may cause cardiac disturbances. Although present in these reagents in very small quantities, appropriate health and safety interventions are recommended. Appropriate Personal Protective Equipment should be worn at all times whilst carrying out tests. Accidental exposure to the skin should be treated by washing with mild detergent and water for 15-20 minutes. Eyes should be washed using water or normal saline for the same duration.

Signed for Quality Control: 

Date: 25 Aug 17

Author: Rafael Carlisle
Technical Reviewer: Paul Booth
QA Reviewer: Jahira Jangaria

SER-DOC-32 v1.0
Authorised Date: 06 Jul 17
Next Review Date: 06 Jul 18

CERTIFICATE OF ANALYSIS

Signed for Certificate Check: 

Date: 25 Aug 18