

Certificate of Analysis

Product name: **Frozen Cell Recovery Solution (FCRS)**

Product code(s): 04-161 – 04-162

Batch Numbers: 41160101

Date Manufactured: 02 Feb 2016

Expiry Date: 01 Feb 2018

Appearance: Clear Liquid

Quality Control:

	Acceptable Ranges	Batch result
pH	6.40-6.80	6.45
Conductivity($\mu\text{s}/\text{cm}$)	44000-52000	48400
Osmolarity(mOsm/kg)	1080-1120	1097

Preservatives used: EDTA (anti-fungal)

Sterility: Terminally filtered with 0.2 μm filter in a grade C clean room and tested for sterility.

Storage: 2-8 $^{\circ}\text{C}$

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

This reagent contains a high concentration of sodium chloride but presents only a very minor risk to staff handling the solution.

Signed for Quality Assurance:

Date: 03 Feb 2016

