

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** LISS (Low Ionic Strength Saline)

**Product code(s):** 04-016 & 04-123

**Batch Number(s):** 03170917

**Date Manufactured:** 25 Sep 2017

**Expiry Date:** 24 Sep 2018

**Appearance:** Clear liquid

**Quality Control:**

	Specification	This batch result
pH	6.50 - 6.90	6.68
Conductivity( $\mu\text{s}/\text{cm}$ )	3500 - 3800	3780
Osmolality(mOsm/kg)	285 - 295	288

**Preservatives used:** Sodium azide (0.09% [w/v])

**Sterility:** N/A

**Storage:** Room temperature

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Wear appropriate Personal Protective Equipment (PPE) when handling this reagent.

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:**



**Date:**

27 Sep 17

**Signed for Certificate Check:**



**Date:**

28 Sep 17