

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Filtered Reagent Grade Water  
**Product code(s):** 04-036  
**Batch Number(s):** 18171015  
**Date Manufactured:** 22 Oct 2017  
**Expiry Date:** 21 Oct 2018  
**Appearance:** e.g. Clear Liquid  
**Quality Control:**

	Specification	This batch result
pH	N/A	N/A
Conductivity( $\mu\text{s}/\text{cm}$ )	0 - 10	5.49
Osmolality(mOsm/kg)	N/A	N/A

**Preservatives used:** N/A

**Sterility:** Product is filtered through a sterile 0.2 $\mu\text{m}$  filter into gamma irradiated bags in a grade B clean room and tested for sterility.

**Storage:** Room Temperature

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

No Health and Safety concerns with this product, however:

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 06 Nov 17

Signed for Certificate Check: 

Date: 06 Nov 17