

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: TBIS - Tris Buffered Isotonic Saline

Product code(s): 04-142

Batch Number(s): 37171007

Date Manufactured: 09Oct17

Expiry Date: 08Oct18

Appearance: Clear Liquid

Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | 7.3-7.6 | 7.37 |
| Conductivity($\mu\text{s}/\text{cm}$) | 14000-16500 | 15080 |
| Osmolality(mOsm/kg) | 285-300 | 294 |

Preservatives used: None

Sterility: Terminally filtered with a 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 26 OCT 17

Signed for Certificate Check:



Date: 26 OCT 17