

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Phosphate Buffered Saline Concentrate (PBSC)
Product code(s): 04-008 - 04-015
Batch Number(s): 04171009
Date Manufactured: 22 Oct 17
Expiry Date: 21 Oct 21
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.95-7.05 (1:25)	7.02 (1:25)
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295 (1:25)	295 (1:25)

Preservatives used: Sulfamethoxazole, Trimethoprim, Sodium azide (0.01% w/v)

Sterility: N/A

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 13 NOV 17

Signed for Certificate Check:



Date: 13 NOV 17