

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: PBS Azide free pH 7.3

Product code(s): 04-047

Batch Number(s): 25171007

Date Manufactured: 20 Oct 17

Expiry Date: 19 Oct 19

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	7.25 - 7.35	7.28
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	286

Preservatives used: N/A

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.


Storage: Ambient

Country of Origin: United Kingdom


Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 14 NOV 17

Signed for Certificate Check: 

Date: 14 NOV 17