

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Modified Alsever's Solution

Product code(s): 04-181, 04-182

Batch Number(s): 32171116

Date Manufactured: 06 Nov 17

Expiry Date: 05 Nov 18

Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.80-7.30	7.16
Conductivity($\mu\text{s}/\text{cm}$)	11500-13000	11747
Osmolality(mOsm/kg)	335-360	345

Preservatives used: Chloramphenicol and Neomycin Sulphate

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: 2-8oC

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date:

16 NOV 17

Signed for Certificate Check:



Date:

16 NOV 17