

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Heart Valve Cryopreservation Medium
Product code(s): 04-311
Batch Number(s): 71171106
Date Manufactured: 06Nov17
Expiry Date: 05Nov18
Appearance: Red translucent liquid

Quality Control:

	Specification	This batch result
pH	7.00-7.50	7.28
Conductivity($\mu\text{s}/\text{cm}$)	8500-12000	9550
Osmolality(mOsm/kg)	2600-3000	2927

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in grade B clean room and tested for sterility


Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 20 Nov 17

Signed for Certificate Check: 

Date: 20 Nov 17