

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Bristol Heart Valve Cryopreservation Medium
Product code(s): 04-342
Batch Number(s): 81171103
Date Manufactured: 10 Nov 17
Expiry Date: 09 Nov 18
Appearance: Red translucent liquid
Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | 7.50-8.50 | 7.59 |
| Conductivity($\mu\text{s}/\text{cm}$) | 1050-1450 | 1149 |
| Osmolality(mOsm/kg) | 1900-2200 | 2017 |

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 20 Nov 17

Signed for Certificate Check:



Date: 21 Nov 17