

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Modified Alsever's Solution

**Product code(s):** 04-182

**Batch Number(s):** 32171117

**Date Manufactured:** 20 Nov 17

**Expiry Date:** 19 Nov 18

**Appearance:** Clear Liquid

### Quality Control:

	Specification	This batch result
pH	6.80-7.30	6.88
Conductivity( $\mu\text{s}/\text{cm}$ )	11500-13000	12537
Osmolality(mOsm/kg)	335-360	346

**Preservatives used:** Chloramphenicol and Neomycin Sulphate

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in a grade B clean room and tested for sterility.

**Storage:** 2-8oC

**Country of Origin:** United Kingdom

### Health and Safety:

#### Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 08 Dec 17

Signed for Certificate Check:



Date: 12 Dec 17