

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Reagent grade water (Non Sterile)
Product code(s): 04-035
Batch Number(s): 17171281
Date Manufactured: 21 Dec 2017
Expiry Date: 20 Dec 2018
Appearance: Clear Liquid

Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | 6.95-7.05 | N/A |
| Conductivity($\mu\text{s}/\text{cm}$) | <10 | 4.43 |
| Osmolality(mOsm/kg) | 285-295 | N/A |

e.g. Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: Reagent water

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: 

Date: 21-Dec-17

Signed for Certificate Check: 

Date: 21-Dec-17