

## CERTIFICATE OF ANALYSIS

# Certificate of Analysis

**Product name:** Heart Valve Cryopreservation Medium  
**Product code(s):** 04-311  
**Batch Number(s):** 71171207  
**Date Manufactured:** 18 Dec 17  
**Expiry Date:** 17 Dec 18  
**Appearance:** Red translucent liquid

**Quality Control:**

	Specification	This batch result
pH	7.00-7.50	7.33
Conductivity( $\mu\text{s}/\text{cm}$ )	8500-12000	9495
Osmolality(mOsm/kg)	2600-3000	2911

**Preservatives used:** None

**Sterility:** Terminally filtered with 0.2 $\mu\text{m}$  filter in grade B clean room and tested for sterility

**Storage:** Ambient

**Country of Origin:** United Kingdom

**Health and Safety:**

**Precautionary Notes**

Please download the Safety Data Sheet at <http://sourcebioscience.com>

**Signed for Quality Control:** 

**Date:** 04 Jan 18

**Signed for Certificate Check:** 

**Date:** 05 Jan 18