

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: e.g. PBSS (Phosphate Buffered Serology Saline)

Product code(s): 04-004
Batch Number(s): 01171220
Date Manufactured: 22 Dec 2017
Expiry Date: 21 Dec 2019
Appearance: Clear Liquid

Quality Control:

	Specification	This batch result
pH	6.95-7.05	6.98
Conductivity($\mu\text{s}/\text{cm}$)	N/A	N/A
Osmolality(mOsm/kg)	285-295	288

e.g. Meets the requirements set out in 'Guidelines for the Blood Transfusion Services in the United Kingdom' 7th Edition; as set out in chapter 12, 12.2.4, Miscellaneous reagents.

Preservatives used: Sodium azide (0.01% [w/v])

Sterility: Product is filtered through a 0.2 μm filter into gamma irradiated bags

Storage: Room temperature

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control: *D. Crossley*

Date: *08-Jan-2018*

Signed for Certificate Check: *Am Green*

Date: *08-JAN-2018*