

CERTIFICATE OF ANALYSIS

Certificate of Analysis

Product name: Amnicryo Medium
Product code(s): 04-326
Batch Number(s): 72180110
Date Manufactured: 03 Jan 18
Expiry Date: 02 Jan 19
Appearance: Red Translucent Liquid
Quality Control:

| | Specification | This batch result |
|---|---------------|-------------------|
| pH | 7.00-7.50 | 7.28 |
| Conductivity($\mu\text{s}/\text{cm}$) | 1250-1600 | 1539 |
| Osmolality(mOsm/kg) | 1680-1850 | 1751 |

Preservatives used: None

Sterility: Terminally filtered with 0.2 μm filter in a grade B clean room and tested for sterility.

Storage: Ambient

Country of Origin: United Kingdom

Health and Safety:

Precautionary Notes

Please download the Safety Data Sheet at <http://sourcebioscience.com>

Signed for Quality Control:



Date: 11 Jan 18

Signed for Certificate Check:



Date: 11 Jan 18